

## Sample preparation

### (1) Chemical

Lysis buffer (10ml)

Urea (9.5 M)	4.8 g
NP-40 (2%)	0.2 ml
Ampholine (2%) (pH 3.5 – 10)	0.2 ml
$\beta$ -mercaptoethanol(5%)	0.5 ml
PVP-40	0.5 g
Water	9.10 ml

Note. This buffer may be stored frozen at  $-20^{\circ}\text{C}$ . Do not continuously freeze and thaw it. Use an aliquot once and discard the remainder

### (2) Method

- Collect sample and homogenize in lysis buffer.
- Centrifuge at 15,000 g for 10 min
- Take the supernatant and centrifuge at 15,000 g for 10 min, again
- Take the supernatant and use for 2D-PAGE